




LANG[®]
Desde 1928




Cardiovascular





FINOCHIETTO

	A	B	C	D	CÓDIGO
NEO NATAL	10 mm	15 mm	75 mm	60 mm	30-006-01
BABY	21 mm	30 mm	30 mm	85 mm	30-006-02
INFANTIL	30 mm	55 mm	55 mm	115 mm	30-006-03
ADULTO	45 mm	65 mm	65 mm	185 mm	30-006-04



FINOCHIETTO ALUMINIUM

FINOCHIETTO BABY

	A	B	C	D	CÓDIGO
NEO NATAL	10 mm	15 mm	75 mm	60 mm	30-006-01A
BABY	25 mm	30 mm	100 mm	85 mm	30-006-02A
INFANTIL	30 mm	55 mm	160 mm	115 mm	30-006-03A
ADULTO	45 mm	65 mm	175 mm	185 mm	30-006-04A



COOLEY ATRIAL

ESQUERDO
30-003-02




COOLEY ATRIAL


DIREITO
30-003-01



HENLEY



14 CM - 14-283-14
17 CM - 14-283-17



14 CM - 14-282-14
17 CM - 14-282-17



16 x 19 mm
16 x 25 mm
16 x 32 mm



N° 1 - 12X200mm
14-300-12

N° 2 - 14X200mm
14-300-14


N° 3 - 16X200mm
14-300-16

N° 4 - 18X200mm
14-300-18


KIRKLIN



DINKHUYSEN
P/AORTA




A	B	CÓDIGO
24 mm	21 mm	30-008-24
28 mm	21 mm	30-008-28
37 mm	27 mm	30-008-37



14 CM - 20-036-14
18 CM - 20-036-18

BAILEY




VALVULA MITRAL
30-008-01



16 CM - 14-240-16
18 CM - 14-240-18

GELPI



WIDIA

17 CM - 20-024-17W
21 CM - 20-024-21W
23 CM - 20-024-23W

AGUDO
14 CM - 14-570-14
16 CM - 14-570-16
18 CM - 14-570-18

ROMBO
14 CM - 14-571-14
16 CM - 14-571-16
18 CM - 14-571-18



WEITLANER



ALICATE
CORTE LATERAL
DUPLA FORÇA



19 CM - 20-102-19
22 CM - 20-102-22

GLUCK



22 CM - 20-120-22

STANDART



350G - 20-380-35
500G - 20-380-50
750G - 20-380-75
1000G - 20-380-10

UNIVERSAL



25 CM
30-660-25

FACA LEBSCHÉ

CARDIOVASCULAR



07 CM - 30-013-07
08 CM - 30-013-08
10 CM - 30-013-10
12 CM - 30-013-12

**BULLDOG
DEBAKEY RETO**



07 CM - 30-014-07
08 CM - 30-014-08
10 CM - 30-014-10
12 CM - 30-014-12

**BULLDOG
DEBAKEY CURVO**

3,5 CM - 30-015-35
5,0 CM - 30-015-50



**BULLDOG
DIEFFENBACH
RETO**

3,5 CM - 30-016-35
5,0 CM - 30-016-50



**BULLDOG
DIEFFENBACH
CURVO**



08mm
CURVO
4,5 CM - 30-018-08
RETO
5,0 CM - 30-017-08

**BULLDOG
DIETRICH**



12mm
CURVO
4,7 CM - 30-018-12
RETO
5,2 CM - 30-017-12

**BULLDOG
DIETRICH**



**PINÇA APLICADORA
DE BULLDOGS**

23 CM
30-110-23



1,2 X 0,6 mm
30-690-12

1,8 X 1,0 mm
30-690-18

2,5 X 1,5 mm
30-690-25

DENNIS



30-692-18



24 CM - 30-702-24

TITANIUM
30-702-20T

DINKHUYSEN
P/ MAMARIA



DEBAKEY



1,5 mm

15 CM - 30-715-15
20 CM - 30-720-15
24 CM - 30-724-15

TITANIUM
1,5 mm

15 CM - 30-715-15T
20 CM - 30-720-15T
24 CM - 30-724-15T

1,8 mm

15 CM - 30-715-18T
20 CM - 30-720-18T
24 CM - 30-724-18T



DEBAKEY



2,0 mm

15 CM - 30-715-20
20 CM - 30-720-20
24 CM - 30-724-20
30 CM - 30-730-20

TITANIUM

2,0 mm

24 CM - 30-724-20T



DEBAKEY




2,5 mm

20 CM - 30-720-25
24 CM - 30-724-25
30 CM - 30-730-25

TITANIUM


2,5 mm

20 CM - 30-720-25T





JAVID
18 CM




4 mm - 30-600-04
6 mm - 30-600-06
8 mm - 30-600-08



GLOVER
16 CM

CURVO
30-030-16

RETO
30-031-16



DEBAKEY
17 CM

30-025-17



DEBAKEY
16 CM

30-025-16



DEBAKEY
16 CM

30-026-16



DEBAKEY
13 CM



45 mm
30-020-12



40 mm - 45°
30-023-11



40 mm - 90°
30-024-10



DeBakey com serrilhas
atraumaticas



DEBAKEY
13 CM



45 mm
30-021-11



DEBAKEY
13 CM




45 mm
30-022-11



DEBAKEY
DERRA



13 CM - 30-060-13
16 CM - 30-060-16
17 CM - 30-060-17
18 CM - 30-060-18



DEBAKEY



17 CM
30-066-17



DEBAKEY
GREGORY



14 CM
30-077-14



14 CM
30-079-14



DEBAKEY



19 CM - 30-084-19
20 CM - 30-084-20



DEBAKEY



20 CM
30-090-20



DEBAKEY
LELAND-JONES



18 CM
30-096-18




19 CM
30-096-19



20 CM
30-096-20






POTTS



POTTS



55 mm
30-112-19

POTTS DEBAKEY
19 CM



55 mm
30-112-20

POTTS DEBAKEY
20 CM



DeBakey com serrilhas
atraumaticas



DEBAKEY

18 CM - 30-180-18
20 CM - 30-180-20
27 CM - 30-180-27
31 CM - 30-180-31



DEBAKEY



16 CM - 30-192-16
17 CM - 30-192-17




DEBAKEY
LAMBERT-KAY



22 CM
30-216-22



21 CM
30-216-21




DeBakey com serrilhas
atraumaticas



18 CM
30-222-18

DEBAKEY
BECK



25 CM - 30-360-25
26 CM - 30-360-26


DEBAKEY
MATOX



DeBakey com serrilhas
atraumaticas



DEBAKEY
60°



DeBakey com serrilhas
atraumaticas



DeBakey com serrilhas
atraumaticas



5 CM - 3 mm
30-010-03

AGULHA C/PONTA
P/HEPARINA



8 mm



10 mm



DESCOLADOR
DE PLACA



DESCOLADOR
DE PLACA



DESCOLADOR
FINO



19 CM
DILATADOR
P/CORONÁRIA





25 CM TOTAL
TUBO 18 CM

30°
30-150-30

CÂNULA
DINKHUYSEN
25 CM



25 CM TOTAL
TUBO 18 CM

90°
30-150-90

CÂNULA
DINKHUYSEN
25 CM



30-666-20

GANCHO
P/SIMPATECTOMIA



30-668-00

PASSA FIO
DE RUMEL



PINÇA
LONGUETTE



23 CM - 13-462-23

DELICADA
23 CM - 13-464-23

CARDIOVASCULAR





34 CM
30-010-00

PINÇA
P/VÁLVULA MITRAL



9 CM
30-770-09

PENTEAN
C/CREMALHEIRA




9 CM
30-772-09

PENTEAN
S/CREMALHEIRA



30-782-18



30-780-18
18 CM
REYNOLD P/TUBOS
SERRILHADA



SACA TROMBO



SEMB CURVA



PORTA AGULHA
PI FIO DE AÇO C/WÍDIA



PORTA AGULHA
MAYO HEGAR
C/WÍDIA



PORTA AGULHA
DEBAKEY C/WÍDIA

- 12 CM - 17-012-12W
- 14 CM - 17-012-14W
- 16 CM - 17-012-16W
- 18 CM - 17-012-18W
- 20 CM - 17-012-20W
- 23 CM - 17-012-23W
- 25 CM - 17-012-25W
- 30 CM - 17-012-30W

- 18 CM - 17-021-18W
- 20 CM - 17-021-20W
- 23 CM - 17-021-23W
- 25 CM - 17-021-25W
- 30 CM - 17-021-30W

- 18 CM
- 17-120-18W
- 20CM
- 17-120-20W



1mm

14 CM - 17-181-14W
16 CM - 17-181-16W
18 CM - 17-181-18W
20 CM - 17-181-20W
25 CM - 17-181-25W

RYDER C/WÍDIA



2mm

14 CM - 17-182-14W
16 CM - 17-182-16W
18 CM - 17-182-18W
20 CM - 17-182-20W
25 CM - 17-182-25W

RYDER C/WÍDIA



18 CM - 17-016-18W
20 CM - 17-016-20W
23 CM - 17-016-23W
25 CM - 17-016-25W
30 CM - 17-016-30W

VASCULAR C/WÍDIA



TITANIUM
18 CM - 17-030-18T
21 CM - 17-030-21T
23 CM - 17-030-23T


BARRAQUER
COM TRAVA C/WÍDIA



TITANIUM
18CM - 17-032-18T
21CM - 17-032-21T
23CM - 17-032-23T

BARRAQUER
SEM TRAVA C/WÍDIA

CARDIOVASCULAR



14 CM - 17-060-14W
18 CM - 17-060-18W
20 CM - 17-060-20W

CASTROVIEJO
COM TRAVA C/WÍDIA



25 CM
17-072-25W

FINOCHIETTO C/WÍDIA



WÍDIA
18 CM - 11-038-18W

BOYD CURVA EM S
VASCULAR



DIETRICH 18 CM



25°
11-420-25

WÍDEA
11-420-25W



45°
11-420-45

WÍDEA
11-420-45W



60°
11-420-60

WÍDEA
11-420-60W



90°
11-420-90

WÍDEA
11-420-90W



125°
11-420-12



POTTS SMITH



25°
18 CM - 11-360-25

WIDEA
18 CM - 11-360-25W



45°
18 CM - 11-360-45

WIDEA
18 CM - 11-360-45W



60°
18 CM - 11-360-60

WIDEA
18 CM - 11-360-60W



METZENBAUM
VASCULAR CURVA


WIDEA

18 CM - 11-306-18W

20 CM - 11-306-20W

23 CM - 11-306-23W

25 CM - 11-306-25W



20 CM - 11-324-20S

TESOURA P/CORTAR
TUBOS SERRILHADA